
HENRY SCHEIN®

900-6368

FINGER PULSE OXIMETER OPERATOR'S MANUAL

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www.henryscheinbrand.com

General Description

Oxygen Saturation is a percentage of Oxyhaemoglobin (HbO₂) capacity, compounded with oxygen, by all combinative haemoglobin (Hb) capacity in blood. In other words, it is the consistency of Oxyhaemoglobin in blood. It is a very important parameter for the Respiratory Circulation System. Many respiratory diseases can result in oxygen saturation being lowered in human blood. Additionally, the following factors can reduce oxygen saturation: Automatic regulation of organ dysfunction caused by Anaesthesia, Intensive Postoperative Trauma, Injuries caused by some medical examinations. These situations might result in light-headedness, asthenia, and vomiting. Therefore, it is very important to know the oxygen saturation of a patient so that doctors can find problems in a timely manner.

The Finger Pulse Oximeter is compact, has low power consumption, convenient operation and portability. It is only necessary for a patient to put one of his fingers into the fingertip photoelectric sensor for diagnosis, and the display screen will show oxygen saturation. It has been proven in clinical experiments that it has high precision and reliability.

Measurement principle

The principle of the Oximeter is as follows: A mathematical formula is established making use of Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive haemoglobin (RHb) and Oxyhaemoglobin (HbO₂) in glow and near-infrared zones.

Operative principle of the instrument: Photoelectric Oxyhemoglobin Inspection Technology is adopted in accordance with Capacity Pulse Scanning and Recording Technology, so that two beams of different wavelengths of light (660nm glow and 940nm near infrared light) can be focussed onto a human nail tip through a clamping sensor. A measured signal obtained by a photosensitive element will be shown on the Oximeter's display through

process in electronic circuits and microprocessor.

Diagram of Operation Principle

1. Red and Infra-red ray Emission Tube

2. Red and Infra- red ray Receipt Tube

Precautions for use

1 Before use, carefully read the manual.

2 Do not use the Finger Pulse Oximeter in an MRI or CT environment.

3 Do not use the Finger Pulse Oximeter in situations where alarms are required.

4 **Explosion hazard:** Do not use the Finger Pulse Oximeter in an explosive atmosphere.

5 The Finger Pulse Oximeter is intended only as an adjunct in patient assessment. It must be used in conjunction with other methods of assessing clinical signs and symptoms.

6 Check the Finger Pulse Oximeter sensor application site frequently to determine the positioning of the sensor and circulation and skin sensitivity of the patient.

7 Do not stretch the adhesive tape while applying the Finger Pulse Oximeter sensor. This may cause inaccurate readings or skin blisters.

8 The Finger Pulse Oximeter has no SpO2 alarms; it is not for continuous monitoring, as indicated by the symbol.

9 Prolonged use or the patient's condition may require changing the sensor site periodically. Change sensor site and check skin integrity, circulatory status, and correct alignment at least every 4 hours.

10 Inaccurate measurements may be caused by autoclaving, ethylene oxide sterilizing, or immersing the sensors in liquid.

11 Significant levels of dysfunctional haemoglobins (such as carbonxy- haemoglobin or methaemoglobin).

12 Intravascular dyes such as indocyanine green or methylene blue

13 SpO₂ measurements may be adversely affected in the presence of high ambient light. Shield the sensor area (with a surgical towel, or direct sunlight, for example) if necessary.

14 Excessive patient movement

15 Venous pulsations

16 Placement of a sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line.

17 The patient has hypotension, severe vasoconstriction, severe anaemia, or hypothermia

18 The patient is in cardiac arrest or in shock.

19 Fingernail polish or false fingernails may cause inaccurate SpO₂ readings.

Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.

Product Properties

1 Operation of the product is simple and convenient.

2 The product is compact, lightweight and fully portable..

3 Low Power consumption.

4 Low voltage warning will be indicated when battery voltage is too low to allow normal operation of the Oximeter.

5 The product will automatically be powered off when inactive for more than 8 seconds.

Product Operation Scope






FINGER PULSE OXIMETER is a portable non-invasive, spot-check, oxygen saturation of arterial haemoglobin (SpO₂) and pulse rate of adult and paediatric patients at home, and hospital (including clinical use in internist/surgery, Anaesthesia, intensive care and etc). It is not for continuously monitoring.

The FINGER PULSE OXIMETER requires no routine calibration or maintenance other than replacement of batteries.

Operation Instructions

- 1 Install two AAA batteries correctly.
- 2 Nip the clamp as in diagram.
- 3 Fully insert one finger into rubber hole of the Oximeter.
- 4 Press the switch button once on the front panel.
- 5 Your finger and body should be steady during measuring.
- 6 Read corresponding data from display screen.

After turning on the Oximeter, each time you press the power switch, the Oximeter will switch to another display mode. There are 6 display modes shown as follows:


When you press and hold the power switch for more than one second, the brightness of the Oximeter will be change there are 10 levels on brightness; the default is level four.

NOTE: Please use medical alcohol to clean the rubber touching the finger on the inside of Oximeter, and clean the test finger using alcohol before and after each use. (The rubber inside of the Oximeter is made from medical rubber, which has no toxin and no harm to the skin).

The nail surface needs to be facing upward when inside the oximeter.

Brief Description of Front Panel

The PR Bar graph displays corresponding with the patient's pulse beat. The height of the bar graph shows the patient's pulse strength.

Product Accessories


1. One lanyard
2. Two batteries
3. One user manual

Battery Installation

1. To open cover - press the button and push in the direction of the arrow.
2. Put the two AAA batteries into battery cassette with correct polarities.
3. To close cover - push down and slide in the opposite direction of the arrow.

Notes:

- Battery polarities must be correctly installed. Otherwise, damage might be caused to device.
- Please put in or remove batteries in carefully, or damage to the device bracket may occur.
- Please remove batteries if the Oximeter will not be used for a long time.

Lanyard Installation

1. Insert thinner end of the strap through the loop.
2. Insert thicker end of the strap through the threaded end before pulling it tightly.

Maintenance and Storage

1. Replace the batteries when low voltage lamp is lit.
2. Clean surface of the fingertip Oximeter before each use/
3. Remove batteries if the Oximeter will not be operated for any length of time..
4. Store the product in a place with ambient temperature $-20^{\circ}55^{\circ}$ and humidity $<93\%$ (no condensation).
5. It is recommended that the product should always be kept in a dry environment. A damp ambient will affect the use or damage the product.
6. Please follow local government guidelines to dispose of used batteries.

Calibration

1. The functional tester cannot be used to assess the accuracy of the Oximeter.
2. Index 2 that made by Bio-Tek company is a function tester. Set Tech to 1, R curve to 2, and then a user can use this particular calibration curve to measure the Oximeter.
3. The test methods used to establish the SpO₂ accuracy is clinical testing. The Oximeter used to measure the arterial haemoglobin oxygen saturation levels and these levels are to be compared to the levels determined from arterial blood sampling with a CO-Oximeter.

Declaration

EMC of this product complies with IEC60601-1-2 standard.

The materials which the user can come into contact have no toxicity and no action on tissues comply with ISO10993-1, ISO10993-5 and ISO10993-10.

Detailed descriptions of product functions

1. **Display Type:** OLED display
2. **SpO₂:** Measurement range: 70-100% Accuracy: 70%-100%, $\pm 3\%$; 69% no definition.
3. **Pulse Rate:** Measure range: 30~235 BPM Accuracy: 30~99bpm, ± 2 bpm; 100~235bpm, $\pm 2\%$
Pulse Intensity: Bar graph Indicator
4. **Power Requirements:**
Two AAA alkaline Batteries Power consumption: Less than 30mA
Low power indication: 
Battery Life:
Two AAA 1.5V, 800mAh alkaline batteries could be continuously operated as long as 30 hours.

5. **Dimension:** Length: 58mm Width: 32mm Height: 37mm Weight: 33g (without batteries)
6. **Environment Requirements:** Operation Temperature: 5°-40° Storage Temperature: -20°55°
Ambient Temperature: <80% in operation, no condensation.
<93% in storage, no condensation.
7. **Measurement Performance in Low Perfusion Condition:**
Required the test equipment (BIO-TEK INDEX Pulse Oximeter tester) the pulse wave is available without failure when the simulation pulse wave amplitude is at 0.6%.
8. **Interference Resistance Capacity against Ambient Light:** Device works normally when mixed noise produced by BIO-TEK INDEX Pulse Oximeter tester.

Guidance and manufacture's declaration – electromagnetic emissions- for all EQUIPMENT and SYSTEMS

Guidance and manufacture's declaration – electromagnetic emission

*The Pulse Oximeter is intended for use in the electromagnetic environment specified below.
The customer or the user of the Pulse Oximeter should assure that it is used in such an environment.*

Emission test Compliance Electromagnetic environment – guidance

RF emissions CISPR 11	Group 1	The <i>Pulse Oximeter</i> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The <i>Pulse Oximeter</i> is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Fault Diagnosis

Problem	Possible reason	Solution
SpO ₂ or PR is not displayed	1. Finger is not inserted correctly 2. Patient's SpO₂ value is too low	1. Re-insert the finger 2. Re-try several times, if no problem with the product then go to a hospital for exact diagnosis
SpO ₂ or PR reading is unstable	1. Finger is not inserted correctly 2. Excessive patient movement	1. Re-insert the finger 2. Restrict patient/ oximeter movement
The Oximeter can not be powered on	1. No battery or low power of battery 2. Batteries might be installed incorrectly 3. The Oximeter might be damaged	1. Replace batteries 2. Reinstall batteries 3. Contact with local customer service centre
Display suddenly turns off	1. The oximeter will automatically power off when no signal is detected longer than 8 seconds 2. Battery power is low	1. Normal 2. Replace batteries
"Error3" or "Error4" is displayed on screen	1. Low power 2. Receiving tube being shielded or damaged 3. Mechanical misplace for receive-emission tube displayed together with broken connector. 4. Amp circuit malfunction	1. Change new battery 2. Please contact with local customer service centre 3. Please contact with local customer service centre 4. Please contact with local customer service centre
"Error7" is displayed on screen	1. Low power 2. Emission tube damaged 3. Current control circuit malfunction	1. Please change battery 2. Please contact with local customer service centre 3. Please contact with local customer service centre

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Henry Schein U.K. Holdings Ltd.