Natural ELEGANCE

UNIVERSAL BOND

Total Etch Light-Cure Adhesive

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Natural Elegance Universal Bond is indicated for bonding to both enamel and dentin after acid etching, in the restoration of teeth.

DESCRIPTION:

Natural Elegance Universal Bond is a light curing adhesive resin system designed to bond to enamel and dentin. Along with Natural Elegance Etching Gel it can provide a quick and easy way to bond composite restorations to the tooth. Natural Elegance Universal Bond provides both priming of adhesion and enamel and dentin.

FORMULACION:

Natural Elegance Universal Bond is a methacrylate based monomer system and photosensitive catalyst in an acetone solvent.

DIRECTIONS FOR USE:

- The tooth should be cleaned by scaling and prophylaxis with flour of pumice. Isolate tooth and prepare cavity form according to normal dental procedures.
- Apply an Etching Gel (we recommend Natural Elegance Etching Gel) following manufacturer's instructions that accompany the product. **Do not allow** the tooth surface to dehydrate after etching and before bonding.
- Dispense 1 – 2 drops of Natural Elegance Universal Bond into a mixing well, or open the unit dose container by twisting the cap and dipping the brush into the bottom of the capsule.
- Keep tooth moist and vigorously brush on one coat of bond. Allow acetone to evaporate.
- Apply 2 more coats of bond. Evaporate the solvent gently with oil-free, water-free air. Do not displace the bonding resin.
- Confirm that all surfaces to be bonded are glossy. If not, reapply Natural Elegance Universal Bond and dry to achieve glossy appearance.
- Using a conventional curing lamp, light cure for 20 seconds with a blue, dental curing light with an output between 470-480 nanometers and a minimum of 500 mW/cm² intensity.
- Proceed with the composite restoration technique.

STORAGE:

- The system is designed to be stored between 35°F-50°F (2°C-10°C).
- Allow equilibration to room temperature prior to use.
- Replace caps immediately after dispensing.
- Do not store materials in proximity of eugenol-containing products.
- Do not store the materials in intense light or under wet conditions.

CONSULT WARNINGS BELOW

Natural ELEGANCE

MICROHYBRID COMPOSITE

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Natural Elegance Microhybrid Composite is indicated for direct esthetic restoration of anterior and posterior teeth, it can be used in limited occlusal function.

DESCRIPTION:

Natural Elegance Microhybrid Composite is a light-cure, resin-based, dental restorative material. It is radio-opaque. It comes in shades matched to the VITA® Shade Guide. The material is a high viscosity composite that will resist sticking to instruments, yet will have sufficient flow when applied to the tooth to adapt to the walls of the cavity preparation. It provides both esthetics and functional durability.

FORMULACION:

Approx. 75w% (54 vol%) silanated 0.7 micron Barium glass and 0.01 micron silica. Approx 25wt % Bis-GMA based dimethacrylate system. Approx 1wt% catalysts, inhibitors and pigments.

DIRECTIONS FOR USE:

- Prophy: The tooth should be cleaned by scaling and prophylaxis with flour of pumice.
- Shade selection: Esthetic choices should be performed prior to isolation and/or preparation of the teeth. Do not stare at the shade tabs and teeth for more than 10-15 seconds during shade selection. Use of a VITA Shade Guide is recommended.
- Isolation: Isolate the teeth. Use of a rubber dam is highly recommended.
- Cavity preparation: Follow usual procedures for tooth cavity preparation, leaving no residual material or base from any previous restoration.
- Pulp Protection: Base the preparation as needed. Calcium hydroxide cement may be used for covering small pulpal exposures. Any additional dentinal coverage desired can be achieved by using a glass ionomer lining cement. Eugenol containing cements must be avoided.
- Placement of Matrix: Place the appropriate matrix or crown form to assure appropriate proximal contour of the tooth. Place wedges to produce proper gingival adaptation and teeth separation. Burnish the matrix to achieve proximal contact. The matrix may be placed following the adhesive application step if preferred.
- Bonding: Use a state of the art dentin/enamel bonding system (we recommend Natural Elegance Universal Bond) to bond this composite to the tooth structure. Modified bonding systems can be used to bond this composite to other kinds of materials used in restorative techniques. Cure the bonding resin before placing the composite. Follow the manufacturer's instructions.
- Dispensing the Composite: Follow the directions corresponding to the dispensing system chosen; **Syringe:** Dispense onto a mixing pad the amount of composite needed to complete the restoration(s). Protect the dispensed material from direct light exposure with an opaque cover. Avoid placing under direct intensity of the operatory light during the delivery process. **Singe Dose Capsules:** Insert the capsule into the Natural Elegance Dispenser for Tips. Deliver an increment of the composite directly into the preparation.
- Placement: The composite should be placed into the preparation in increments no greater than 2.5 mm. **It is important to use thinner layers for darker shades.** Place the composite into the cavity using an appropriately designed dental hand instrument (some metal instruments may leave gray marks on the composite) and adapt the composite to the walls of the preparation. Should the material begin to stick to the instrument,

the instrument can be wiped clean with a dry, alcohol gauze. Do not allow any contamination or washing between cured layers.

- Curing: Cure each layer for 30 seconds with a blue dental curing light with an output between 470-480 nanometers and a minimum of 500 mW/cm² intensity. An additional 20 seconds may be used from other directions. Any wedge and/or matrix can be removed after the last layer of composite is cured. Do not allow any contamination with saliva or washing between cured layers.

- Finishing: After curing, remove flash with a sharp instrument. If necessary adjust restoration to proper margins, contours, contacts and occlusion. Finish using a fine diamond, multifluted carbide finishing bur, or finishing disc and strips.
- Polishing: Various polishing pastes are available which can enhance the polish.
- Resealing: After the polishing procedure is completed, microscopic surface defects and marginal bonding that is stressed can be resealed. The composite surface can be cleaned and the surrounding enamel re-etched with phosphoric acid, washed and dried thoroughly. The bonding agent is applied in a thin layer, dried and cured.

STORAGE:

- The system is designed to be stored between 50°F-75°F (10°C-24°C).
- Replace caps immediately after dispensing.
- If stored in a refrigerator, bring to room temperature prior to use.
- Do not store materials in proximity of eugenol-containing products.
- Do not store the materials in the intense light, or under wet conditions.

CONSULT WARNINGS BELOW

Natural ELEGANCE

FLOWABLE COMPOSITE

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Natural Elegance Flowable Composite is indicated for class III and V restorations. It can also be used for repair of crown margins at the discretion of the dentist. In posterior composite restorative techniques this flowable composite may also be used immediately after the bonding adhesive at the base of the restoration. It can be used to modify the margins of (bisacryl) temporary crowns. It can provide solutions to a number of clinical problems in dental restorative technique especially where flow can improve adaptation.

DESCRIPTION:

Natural Elegance Flowable Composite is a light-cure, resin based, dental restorative material. It is radio-opaque. It comes in shades matched to the VITA® Shade Guide. Its low viscosity and precise syringe delivery make it flowable, easy to place and finish, and adaptable to a variety of clinical situations.

FORMULACION:

Natural Elegance Flowable Composite is a 41% filled (by volume) composite based on a resin system of BIS-GMA and triethylene glycol dimethacrylate. The filler has particle sizes ranging from submicron to 8 microns. Contains 2% Sodium Fluoride (by weight). There is about 1wt% catalysts, inhibitors and pigments.

DIRECTIONS FOR USE:

- Prophy: The tooth should be cleaned by scaling and prophylaxis with flour of pumice.
- Shade selection: Esthetic choices should be performed prior to isolation and/or preparation of the teeth. Do not stare at the shade tabs and teeth for more than 10-15 seconds during shade selection. Use of a VITA Shade Guide is recommended.
- Isolation: Isolate the teeth. Use of a rubber dam is highly recommended.
- Cavity preparation: Follow usual procedures for tooth cavity preparation, leaving no residual material or base from any previous restoration.
- Pulp Protection: Base the preparation as needed. Calcium hydroxide cement may be used for covering small pulpal exposures. Any additional dentinal coverage desired can be achieved by using a glass ionomer lining cement. Eugenol containing cements must be avoided.
- Bonding: Use a state of the art dentin/enamel bonding system (we recommend Natural Elegance Universal Bond) to bond this composite to the tooth structure. Modified bonding systems can be used to bond this composite to other kinds of materials used in restorative techniques. Cure the bonding resin before placing the composite. Follow the manufacturer's instructions.
- Dispensing the Composite: Attach a syringe tip to the Natural Elegance Flowable Composite syringe. The composite can be applied directly to the tooth surface by exerting gentle pressure on the syringe plunger. Stop pressure on plunger before the desired amount of material is dispensed, since the flow of material from the syringe tip continues for a short time after the pressure is released. Slight back pressure on the plunger can help to decrease the continued flow of the resin out of the syringe.
- Placement and Curing: A layer of flowable composite should not exceed 2 to 2.5 mm. Light cure each layer for 30 seconds with a blue, dental curing light with a wavelength between 470-480 nanometers and a minimum of 500 mW/cm² intensity.
- Finishing: After curing, remove flash with a sharp instrument. Natural Elegance Flowable Composite provides excellent smoothness and may simplify finishing. If needed, finish using a fine diamond, multifluted carbide finishing bur, or finishing disc and strips.
- Polishing: Various polishing materials are available that give excellent results.

STORAGE:

- The system is designed to be stored between 50°F-75°F (10°C-24°C).
- Replace caps immediately after dispensing.
- If stored in a refrigerator, bring to room temperature prior to use.
- Do not store material in proximity to eugenol-containing materials.
- Do not expose to elevated temperatures or intense light.

CONSULT WARNINGS BELOW

Natural ELEGANCE

SEALANT

Light-Cure Opaque Pit and Fissure Sealant

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Natural Elegance Sealant is indicated for the covering and sealing of pits and fissures in the occlusal surfaces of posterior teeth to aid in the resistance to caries in those surfaces.

DESCRIPTION:

Natural Elegance Sealant is a light-cured pit and fissure sealant. It is a very low viscosity resin system with a white opaque color to help identify its retention at recall appointments. The direct delivery

syringe system provides a quick and easy application method. The resin cures quickly to a hard and retentive sealant.

FORMULACION:

Natural Elegance Sealant is a dilute BisGMA based resin system with a white opaquer, a microfine silica particle and a photosensitive catalyst system.

DIRECTIONS FOR USE:

- Clean the teeth with a non-fluoride, oil-free prophy paste. Wash with an oil-free air/water spray.
- Isolate teeth with cotton rolls or rubber dam.
- Blow dry completely with an oil/moisture-free air source.
- Dispense the etchant (we recommend Natural Elegance Etchant) into a mixing well and apply the etchant to the occlusal grooves of the teeth with a brush. Rubsc bottle after dispensing. Agitation of the etchant with the brush bristles may help to increase flow into the fissures. Allow etchant to remain for 30 seconds.
- Wash area completely for 20 seconds with adequate suction to remove etchant and water. **Do not allow saliva to contaminate etched enamel surfaces.**
- Air dry thoroughly using oil/moisture-free air.
- Etched surfaces of teeth should have a dull, white, frosted appearance when etched and dried. If the enamel is not chalky appearing, repeat the process allowing etchant to remain on the teeth for additional 30 second increments of time and wash and dry. Several repetitions of this etching procedure may be needed especially on primary teeth and teeth high in fluoride content.
- Apply a syringe tip to the end of the Natural Elegance sealant syringe. (After completing treatment on a patient, the syringe tip should be removed and properly disposed. Replace the original cap on the syringe after patient treatment.)
- Apply Natural Elegance sealant directly onto each etched, dry occlusal surface taking care to avoid air entrapment. After dispensing the sealant on the teeth, slight back pressure on the plunger will prevent continued flow of the resin out of the syringe. Agitation of the unset sealant with a brush may help adaptation into the pits.
- Expose all coated areas to a dental light curing source for 30 seconds each. Use a blue curing light (bandwidth between 470-480 nanometers with a minimum of 500 mW/cm² intensity.) Hold the end of the light guide as close to the tooth as possible without touching the resin.
- After curing, wipe the occlusal surfaces with a cotton roll to remove the oxygen inhibited layer at the surface of the resin, or rinse with copious amounts of water with adequate suction.
- Check for complete coverage of the fissures by using an explorer. Incomplete coverage requires re-etching and reapplication of sealant. Adjust occlusal surface to avoid occlusal interferences.

STORAGE:

- The material is designed to be stored between 35°F-50°F (2°C-10°C). Allow equilibration to room temperature prior to use.
- Replace caps immediately after dispensing.
- Do not store materials in proximity of eugenol-containing products.
- Do not store the materials in intense light or under wet conditions.

WARNINGS:

- Test the dental curing light periodically to assure sufficient intensity of the light output.
- Dentists and assistants should wear gloves, masks and protective eye wear. Patients also should wear eye protection.
 - If contact with eyes occurs, flush with copious amounts of water.
 - If contact with skin occurs, wash with soap and water.
 - In case of ingestion, seek medical attention immediately.
- Etchant gels are caustic and can cause tissue damage in contact with skin or eyes.
- Uncured monomers in this resin material may cause an allergic reaction or skin irritation in some individuals. Do not use in individuals with allergies to these materials.
- These products contain methacrylate resins. Avoid use of these products on patients with known methacrylate allergies. To reduce the risk of allergic response, avoid exposure to uncured resin.
- Any non-sterilizable item that is handled in the delivery of the dental service should be disinfected by standard dental office hygiene procedures.
- Do not use eugenol-containing materials for pulp protection since they can retard the curing process.
- After applying Natural Elegance Universal Bond, do not allow saliva or water contamination of unset or set material, or of any step in the subsequent restorative procedure. Universal Bond is flammable.
- With Natural Elegance Sealant, if a fluoride treatment is planned at the same appointment, apply the fluoride after bonding the sealant and the sealant is completely set. Do not allow saliva or water contamination of etched tooth surface or unset material.

- For additional information consult the Safety Data Sheet.

HENRY SCHEIN®

Natural ELEGANCE

Distributed by (in US only)
HENRY SCHEIN INC.
135 DURYEA ROAD
Melville, NY 11747 USA

EC REP
HENRY SCHEIN U.K. HOLDINGS LTD.
MEDCARE HOUSE, CENTURION CLOSE
GILLINGHAM BUSINESS PARK
GILLINGHAM, ME8 OSB U.K.

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