Henry Schein Vinylpolysiloxane (VPS Hydro) Impression Material Putty and Monophase for automatic mixing machine (5 :1 ratio) – Instructions for use

Description:

ISO 4823 Type 0 Putty Consistency - ADA 19 Type 1 Very High Viscosity

ISO 4823 Type 2 Medium-Bodied Consistency - ADA 19 Type 1 Medium Viscosity

Indications:

VPS Hydro Putty 380 is indicated for filling the impression tray in the single-stage technique with two materials in combination with VPS Hydro Light Body, Regular Set. VPS Hydro Mono 380 is indicated as single viscosity material for special tray in the monophase technique with one material.

Impression procedure:

The cylinder fits the automatic mixing machines. Refer to the separate instructions supplied by the machine manufacturer.

Carefully remove the cartridge from the package. Remark: The cartridge is heavy. Do not let it fall as it may be damaged and become unusable. – Firmly hold the cartridge with dispensing nozzles up. — Firmly grasp the tab while supporting the cartridge, and remove it from the cylinder moving it upwards. Do not place undue stress on tab. Remark: Insert the cartridge in the mixing machine according to the instructions supplied by the machine manufacturer. – Extrude a small quantity of materials from the dispensing nozzles at the same time (This should be done prior to each use). – Wipe off the extruded material in a vertical direction to avoid cross contamination of the dispensing nozzles. – Install a dynamic mixer on the dispensing nozzles. Remark: Do not use other manufacturer’s dynamic mixers as they are not compatible. – If the dynamic mixer does not fit easily, be sure the internal central hexagonal section of the mixer is aligned correctly with the drive shaft of the mixing machine. – Once the dynamic mixer is properly placed, place the bayonet ring on the dynamic mixer and turn to the right (1/4 turn) until ring is locked in place. – Start the automatic mixing machine up and fill a suitable tray (Putty and Monophase) or an intraoral syringe (Monophase only). – After each use take the cartridge off the machine and stand upright keeping the dynamic mixer inserted.

Disinfection:

Impression may be disinfected or sterilized with a disinfectant for impressions.

Casting:

Store impression at room temperature. Models may be cast immediately after setting.

Compatibility with gypsum:

These products are compatible with all well known gypsum materials (type 3 Gibraltar Stone – Henry Schein; type 4 Gibraltar Die Stone – Henry Schein).

Electroplating:

Impressions may be copper or silverplated.

Storage:

Do not store above 27°C (80°F). Shelf life 3 years when stored at 23°C (73°F).

General notes:

Avoid direct contact with latex gloves as they may inhibit polymerization.

TECHNICAL DATA	VPS HYDRO PUTTY 380 Regular Set	VPS HYDRO PUTTY 380 Fast Set	VPS HYDRO MONOPHASE 380
Mixing Ratio Base : Catalyst	5 : 1	5 : 1	5:1
Mixing Time	Auto Mixing	Auto Mixing	Auto Mixing
Total Working Time (incl. Mixing Time)	2 min.	1 min. 30 secs.	2 min.
Time in mouth (minimum)	3 min. 30 secs.	2 min. 30 secs.	3 min. 30 secs.
Setting Time	5 min. 30 secs.	4 min.	5 min. 30 secs.
Linear Dimensional Change after 24 Hrs.	<-0.20%	<-0.20%	<-0.20%
Maximum Strain in Compression	3%-5%	3%-5%	3%-5%
Recovery from Deformation	>99,5%	>99,5%	>99,5%
Hardness « Shore A »	60±2	60±2	60±2

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Henry Schein U.K. Holdings Ltd.