
CE

 HENRY SCHEIN®

900-3251

ALGINATE MIXER

MEZCLADOR DE ALGINATO
MÉLANGEUR POUR ALGINATE
HS-ALGINATMIXER
MISCELATORE DI ALGINATO

INSTRUCTIONS FOR USE
INSTRUCCIONES DE USO
MANUEL D'UTILISATION
GEBRAUCHSANWEISUNG
GEBRUIKSAANWIJZING

Distributed by:
Distribué par:

Distribuido por:
Vertrieb durch:
Distribuito da:


HENRY SCHEIN INC.
Melville, NY 11747 USA

Henry Schein U.K. Holdings Ltd.
Gillingham ME8 0SB U.K.

Users manual for ALGINATE MIXER

Table of Contents

INTRODUCTION	1
EXCLUSIONS.....	1
ADVISORY WARNINGS to prevent accident or malfunction	1
DETAILS AND ILLUSTRATIONS of components and accessories...	2
USER ADVICE.....	2
ENVIRONMENTAL CONDITIONS for Installation and Use...	3
INSTRUCTIONS.....	4
SET THE MIXING TIME	5
MAINTENANCE	6
TROUBLESHOOTING.....	6
SPECIFICATION.....	7
CHECK LIST OF SERVICE & MAINTENANCE.....	7

INTRODUCTION :

Thank you for purchasing the HS Alginate Mixer. The HS Alginate Mixer is a unique apparatus for mixing dental alginate used by dental practices and orthodontic professionals. Designed with built-in microprocessors and pre-programmed, the HS Alginate Mixer gives the highest performance efficiency to provide a complete bubble-free alginate mix of a fine consistency to obtain the perfect impression.

The HS Alginate Mixer uses centrifugation and by combining vibration and rotation mixes the dental alginate impression material in accordance with the ratio of powder/water and the pre-set recommended mixing time. None of items described in this User Manual apply to any other materials.

Prior to using the equipment, please read the **USER MANUAL** carefully.

Please keep the **USER MANUAL** in a safe place for future reference. Any modifications to subsequent models will not be notified to present users.

The Henry Schein Alginate Mixer is for dental use only; please ensure that the apparatus is operated in accordance with the instructions in the user manual.

Any damage or repairs should be carried out by a fully trained professional. Please contact your local distributor.

EXCLUSIONS :

1. We shall not assume any responsibility for any malfunction, damage, accident or personal injury caused by the improper removal, modification, maintenance or repair not being carried out by Henry Schein or their representative.
2. We shall not assume any responsibility for any malfunction, damage, accident or personal injury caused by the use of product not provided by us.
3. We shall not assume any responsibility for any malfunction, damage, accident or personal injury caused by the improper removal, modification, maintenance or repair with parts not provided by us.
4. We shall not assume any responsibility for any malfunction, damage, accident or personal injury caused by the use not in compliance with the manual.
5. We shall not assume any responsibility for any malfunction, damage, accident or personal injury caused by changes to the ambient environment not in compliance with the power installation and user conditions specified on this manual.
6. We shall not assume any responsibility for any malfunction, damage, accident or personal injury caused by fire disaster, earthquake, flood, or lightning strike.

ADVISORY WARNINGS to prevent accident or malfunction :

Before use please read the "Important warning notices" marked below. All operation should be in compliance with any important notices specified in the manual to avoid any malfunction.

WARNING NOTICES

- ⚠ **DANGER:** indicates the risk of causing death or serious injury if the product is operated incorrectly without due care and attention.
- ⚠ **WARNING:** indicates the risk of causing malfunction or damage to the Henry Schein Alginate Mixer if the product is operated incorrectly without due care and attention.
- ⚠ **CAUTION:** indicates the risk of causing a potential slight injury to the user and other persons, or damage to the product without due care and attention.
- ! **IMPORTANT:** indicates risk of malfunction of machine, imperfect results, poor mixture or other effects.

DETAILS AND ILLUSTRATIONS components and accessories :

1. Main switch
2. Power on/off button
3. Memory selection
4. Start/Stop button
5. Mixing Time display
6. Increase▲/Decrease▼ Timer set
7. Memory address indicator
8. Power cord
9. Voltage adaptor (100/200~240 Vac.)
10. Fuse holder
11. Mixing bowl
12. Machine casing
13. Spatula
14. Mixing stick
15. Water-dosing bottle
16. Basket for miscellaneous articles
17. Retention bolts (immobilisation bolts) in base of unit.

USER ADVICE :

Please adhere fully to the following operational guidelines to ensure smooth, malfunction free operation.

1. Use cold water (10C° ~ 20C°) to extend the setting time. If using a fast-set alginate monitor the setting time carefully.
2. Before operating ensure that the Mixing bowl lid has been fully closed, to avoid any spillage or malfunction.
3. Capacity of alginate impression material for mixing is between 25g and 100g, never exceed 100g.
4. Only operate with a closed mixing cup to prevent the possible spillage while running.
5. Ensure small hole in the centre of Mixing bowl lid is free of debris or residue or this will impair closing.

6. Should any abnormal noise or vibration occur while in mixing, press the START/STOP button immediately and lift the cover to check contents.
7. The Mixing bowl is replaceable. Renew if damage or difficulty in closure is experienced.


Fig.1

8. There are two retention bolts (for immobilisation) at the base of the unit, these should be removed before use and replaced if the unit is moved or sent for repair, (Fig.1) this prevents unnecessary damage during transit.

ENVIRONMENTAL CONDITIONS FOR INSTALLATION AND USE :

1. Conditions for installations and use:

- 1) For indoor use only.
- 2) Operating ambient temperature is between 2°C ~ 40°C (35.6 ~ 104°F).
- 3) Operating altitude is below 2000 metres and atmospheric pressures between 900 ~ 1600 millibars.
- 4) Relative Humidity is lower than 80%.
- 5) Operating area should be free of dust, acids or solvent vapours.
- 6) The required electric voltage: AC110V (+/- 15V); AC220 (+/- 20V).

2. Installation:

- 1) The Henry Schein Alginate Mixer should be positioned on a firm, flat, vibration free surface.
- 2) Never install near to inflammable substances.
- 3) Keep away from water or moisture to avoid electric shock or mechanical problems.
- 4) Do not place anything on the cover or control panel to avoid any damage or spillage.

3. Before/After Transportation:

- 1) Before installing the Henry Schein Alginate Mixer remove the two bolts fixed at the bottom and follow the instruction (Fig.1) to store them safely.
- 2) If the mixer needs to be moved or sent for repair, replace the bolts to immobilize the motor to prevent possible malfunction or damage.

4. Power:

- 1) Check the power switch on the back of the unit is off, and then connect the power cord to the voltage adaptor. (Fig. 2)
- 2) Plug into a standard AC outlet, with co-axial cable of AC 220V 50/60Hz with a 3Amp fuse or AC110V 50/60Hz with a 5Amp fuse. In order to avoid electric shock, the electric power cord shall be earthed.

Fig.2

⚠️ DANGER:

- Check to see the Henry Schein Alginate Mixer's voltage requirement matches your power source. Damage to the Henry Schein Alginate Mixer, electric shock, or fire may occur if voltages are incompatible.
- Keep any potentially dangerous substances a minimum of 30cm (1ft) radius from the Henry Schein Alginate Mixer.
- Ensure that hand and body are dry when touching power socket, main switch or operating the equipment.

INSTRUCTIONS :

Mixing alginate:

- 1) Turn on the main switch on the back of the machine (Fig.3) An LED indicator will light on the left side of the control panel and the Henry Schein Alginate Mixer is ready for operation.
- 2) Press "POWER" on/off button on the control panel to display time indicator showing digital pre-set mixing time (Fig.4)

Fig.3

Fig.4

CAUTION:

- Always keep the equipment in a sleep mode with "the time display is off" when not in use.
- Do not operate empty. The Henry Schein Alginate Mixer is designed to operate when material is in the Mixing bowl. **No impression material, no operation!**

- 3) Measure the appropriate amount of powder using the scoop and instructions supplied with the chosen alginate material.

NOTE: Before dosing, agitate the alginate powder to loosen it well. Scoop it from the storage box with a light, swift movement and then stroke the powder smooth: do not compress the powder in the scoop.

- 4) Using the water dosing bottle or a beaker measure the appropriate amount of cold water at room temperature (10~20°C). Add water to the alginate in Mixing bowl.

IMPORTANT: The setting time of alginate can be controlled by different water temperature. Care should be taken as setting time will be shortened if cold water is not used, particularly in summer.

Fig.5

Fig.6

- 5) Follow to the instruction of the recommended mixing ratio for a better consistency of impression. You can make the mixture thinner (higher viscosity) by adding more or less water respectively (Fig.5).
- 6) Stir the alginate powder and water together in the Mixing bowl for about 2 seconds before placing it into the mixing well (Fig.6). The user can alternatively, shake the cup a few times for preliminary mixing before placing it into the mixer if not using the stirrer provided. Skipping this pre-mix step may result in poor mixing or imperfect consistency of mixture.

Fig.7

Fig.8

- 7) Close the lid of the Mixing bowl tightly. Ensure the small hole in the centre of the lid is not blocked. (Fig.7/Fig.8)

HINT: Scrape the stirrer along the edge of mixing bowl to remove excess material and avoid any residual staining.

- 7) Close the lid of the Mixing bowl tightly. Ensure the small hole in the centre of the lid is not blocked. (Fig.7/Fig.8)

IMPORTANT: The Mixing bowl lid must be closed tightly; otherwise spillage or splash may occur or separated lid may hit against the cover. The Mixing bowl lid is replaceable, please renew if it becomes loose fitting.

Turn

8) Place the Mixing bowl securely into the mixing well. (Fig.9)
NOTES: Follow these instructions swiftly to avoid the early setting of alginate mixture. Remember, longer preparation, leaves shorter setting and working time.

Fig.9

9) Machine will not operate unless the lid is properly closed.
 10) Select the mixing time from the table in this Page for the quantity of alginate used. Press the **START** button to run operation (Fig.10). For safety, the cover of the mixer is locked during operation. The digital timer will count down in 1 second intervals and at completion a “beep” will sound and the cover will unlock.

Fig.10

IMPORTANT: Changing the mixing time will affect the setting and working time of the alginate. Please set the mixing time as short as required to obtain longer setting and working time.

NOTE: The working/setting time increases with a lower (water) temperature. At higher temperatures the time becomes shorter. The water hardness has the same effect; harder the water used for mixing, the shorter the working/setting time.

Fig.11

IMPORTANT: Pressing the button **START/STOP** will stop the mixer immediately. If any abnormal conditions occur during operation, the mixer should be also stopped immediately to conduct an inspection. If the abnormal condition cannot be rectified, please return mixer for repair.

11) When mixing is complete, open the cover and remove the Mixing bowl. Use the spatula to remove the alginate mixture and fill the required impression tray(s) immediately. (Fig.11)

WARNING: Although the mixer stops automatically, do not open the cover until the “beep” has sounded.

SET THE MIXING TIME :

(1) Memory Pre-programmed	Mixing Time (sec.)	(2) Scoop of Alginate	Mixing Time (sec.)
1	8	1 (25g incl. water)	7~10
2	10	2 (35g incl. water)	8~11
3	12	3 (45g incl. water)	9~12

1) The timer has three pre-set memory programs (1).

(Each memory can be altered and re-programmed by the user).

2) Press **MEMORY/SELECT** button to select time required for mixing.

3) To manually adjust the time press the **▲ or ▼** button. The maximum mixing time is 16 seconds. When the mixing is complete, the mixing time will return to its original pre-programmed time.

4) To permanently modify a pre-set mixing time, press the **▲ or ▼** button to the required mixing time. Then press **▲ or ▼** simultaneously. The mixer will sound 2 beeps and mixing time display will flash twice.

5) The new mixing time will now be stored as a pre-set in the memory.

CAUTION: Increase mixing time hastens the setting of the impression material.

*If using 50Hz (with water temperature of 10-20°C) time should be set as shown above (2).

**If using 60Hz (with water temperature of 10-20°C) time should be set as shown above (3).

IMPORTANT:

1. The above recommendation re alginate measurement and mixing time is an average time. Any difference in water temperature or room temperature suggested by the alginate manufacturer will require a different mixing time.

2. Always adhere to the manufacturers instructions when set the mixing time for the quantity used.

MAINTENANCE :

1. Turn off the main switch and disconnect the Henry Schein Alginate Mixer from its power source before any cleaning or maintenance procedures.
2. Do not clean while the mixer is operating or at high temperatures.
3. **Mixing bowl:** After setting the remained residues of impression material inside the Mixing bowl can be cleaned and removed easily. Always remove any blocks or residue from the small hole at the centre of the lid. (Fig.12)


Fig.12

! IMPORTANT:

Cleaning: Residues on the inside of the mixer or mixing bowl or inside of the lid can affect the mixing power. Clean these areas with a dry cloth or water-soaked cloth (not wet). If serious stained, please use a neutral cleaning agent. Ensure any liquid used in cleaning is removed before next use as it may cause problems on subsequent mixing.

⚠ CAUTION:

Never use the solvent or volatile oil for cleaning, this may damage the machine.

TROUBLESHOOTING :

! **IMPORTANT:** Below is a list of hints re problem solving. If problems still persist, please contact Henry Schein and return item for repair.

PROBLEM	CAUSES	SOLUTIONS
Henry Schein Alginate Mixer does not turn on	<ul style="list-style-type: none"> • Power cord is disconnected from power source • The Fuse is blown 	<ul style="list-style-type: none"> • Inspect and connect the power cord to the outlet • Inspect and/or change the fuse
Henry Schein Alginate Mixer does not work	The lid is not closed	Close the lid tightly
Abnormal noise or vibration in running	<ul style="list-style-type: none"> • The screw nuts have not been removed from the bottom of the unit • Start but idlerunning nothing inside the Mixing bowl • The amount of mixture is too much • Mixing bowl lid come off • Damage of gears 	<ul style="list-style-type: none"> • Remove 2 screw nuts from the bottom • Stop the mixer. Make sure to add the mixture into the cup • Limit the mixture to 3 scoops • Stop the mixer. Close Mixing bowl lid well • Send for repair
Time display does not lit on	PCB faulty	Send for repair
Poor mixing	<ul style="list-style-type: none"> • Incorrect ratio of powder/water • Preliminary mixing is not done (for bulk amount mixture over 2 scoops) • Too short of mixing time • Insufficient speed of rotation • Poor quality of alginate or expired 	<ul style="list-style-type: none"> • Check correct ratio of powder / water • Pre-mix the alginate powder in the Mixing bowl with stirrer, or shake a few times, before put into the mixer • Reset the mixing time • Send for repair • Change a new batch alginate
Error code "E1" appear on the time display	The top cover is not properly closed	Close the top cover again
Error code "E2"	Defective control board	Send for repair

SPECIFICATION :

Model:	Henry Schein Alginate Mixer
Power:	AC 110V (+/-15V), AC220V (+/-20V), 220W
Maximum Speed:	3400 r.p.m (60Hz) 2900 r.p.m (50Hz)
Motor:	AC Motor, 2 pole (110V 220W; 220V 200W)
Dimensions:	205 (W) x 250 (D) x 290 (H) mm
Weight:	41.9 lb (19 kg)
Capacity:	25~100g (bowl)
Timer:	1~16 seconds, adjustable
Programs:	3 Memories (pre-set 8, 10 & 12 seconds)
Safety functions:	Automatic door locking and stopping
Operating Conditions:	Temperature: 2°C~40°C (35.6~104°F)
Humidity:	lower than 80%
Altitude:	below 2000 meters

CHECK LIST OF SERVICE & MAINTENANCE :**NOTE:**

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