

Itemcode	Description / Applicable Area	Contents	Classification Medical Device Class / Biocide / Cosmetical / Technical product	Bacteria in accordance to	Candida albicans in accordance to	TbB in accordance to	Membrane covered viruses i.e. HBV, HIV, HCV) in accordance to	Rotavirus	MRSA	H1N1	H1N5	Others
9001288	Eurosept Max Bur	2l	MD IIA	VAH/DGHM, EN1040	VAH/DGHM, EN 1275	VAH/DGHM,	BGA(RKI)/DVV.	—	x	—	—	Adeno Viruses acc, DW
9001289	Eurosept Max Bur	5l	MD IIA	VAH/DGHM, EN1040	VAH/DGHM, EN 1275	VAH/DGHM,	BGA(RKI)/DVV.	—	x	—	—	Adeno Viruses acc. DVV
9001294	Eurosept Max Evac	2l	MD IIA	VAH/DGHM, EN 1040	VAH/DGHM EN 1275	—	BGA(RKI)/DVV and EN 14476	—	x	—	—	—
9001295	Eurosept Max Evac	5l	MD IIA	VAH/DGHM, EN 1040	VAH/DGHM and EN 1275	—	BGA(RKI)/DVV and EN 14476	—	x	—	—	—
9002098	Eurosept Max Impression Liquid	5l	MD IIA	VAH/DGHM	VAH/DGHM	—	BGA(RKI)/DVV.	—	x	—	—	—
9001296	Eurosept Max Impressions powder	700g	MD IIA	VAH/DGHM	VAH/DGHM	VAH/DGHM	BGA(RKI)/DVV.	—	x	x	—	Adeno Viruses acc. DVV
9003331	Eurosept Max Instrument, blue	2l	MD IIA	VAH/DGHM,EN 13727, EN 14561	VAH/DGHM, EN 13624, EN 14562	VAH/DGHM	BGA(RKI)/DVV. EN 14348/ 14563	x	x	x	—	Adeno Viruses acc. DVV
9003332	Eurosept Max Instrument, blue	5l	MD IIA	VAH/DGHM,EN 13727, EN 14561	VAH/DGHM, EN 13624, EN 14562	VAH/DGHM	BGA(RKI)/DVV. EN 14348/ 14563	x	x	x	—	Adeno Viruses acc. DVV
9002600	Eurosept Max Instrument Lemon	2l	MD IIA	VAH/DGHM,EN 13727, EN 14561	VAH/DGHM, EN 13624, EN 14562	VAH/DGHM	BGA(RKI)/DVV. EN 14348/ 14563	—	x	—	—	—
9002601	Eurosept Max Instrument Lemon	5l	MD IIA	VAH/DGHM,EN 13727, EN 14561	VAH/DGHM, EN 13624, EN 14562	VAH/DGHM	BGA(RKI)/DVV. EN 14348/ 14563	—	x	—	—	—
9001291	Eurosept Max Surface	1l	MD IIA	VAH/DGHM, EN 1040, EN 13727	VAH/DGHM, EN 1275, EN 13624	EN 14348 (high bio load)	BGA(RKI)/DVV	x	x	x	x	Adeno Viruses (EN 14476)
9001290	Eurosept Max Surface	5l	MD IIA	VAH/DGHM, EN 1040, EN 13727	VAH/DGHM, EN 1275, EN 13624	EN 14348 (high bio load)	BGA(RKI)/DVV	x	x	x	x	Adeno Viruses (EN 14476)
9002092	Eurosept Max Surface Concentrate	5l	MD IIA	VAH/DGHM,	VAH/DGHM	EN 14348 (high bio load)	BGA(RKI)/DVV	x	x	x	x	Adeno Viruses (EN 14476)
9002089	Eurosept Max Surface Lemon	1l	MD IIA	VAH/DGHM, EN 1040, EN 13727	VAH/DGHM, EN 1275, EN 13624	EN 14348 (high bio load)	BGA(RKI)/DVV	x	x	x	x	Adeno Viruses (EN 14476)
9002090	Eurosept Max Surface Lemon	5l	MD IIA	VAH/DGHM, EN 1040, EN 13727	VAH/DGHM, EN 1275, EN 13624	EN 14348 (high bio load)	BGA(RKI)/DVV	x	x	x	x	Adeno Viruses (EN 14476)
9001293	Eurosept Max Wipes	150	MD IIA	VAH/DGHM, EN 1040, EN 13727	VAH/DGHM, EN 1275, EN 13624	EN 14348 (high bio load)	BGA(RKI)/DVV	x	x	x	x	Adeno Viruses (EN 14476)
9001292	Eurosept Max Wipes	150	MD IIA	VAH/DGHM, EN 1040, EN 13727	VAH/DGHM, EN 1275, EN 13624	EN 14348 (high bio load)	BGA(RKI)/DVV	x	x	x	x	Adeno Viruses (EN 14476)
9002095	Eurosept Max Wipes Lemon XL	200pk	MD IIA	VAH/DGHM, EN 1040, EN 13727	VAH/DGHM, EN 1275, EN 13624	EN 14348 (high bio load)	BGA(RKI)/DVV	x	x	x	x	Adeno Viruses (EN 14476)
9002096	Eurosept Max Wipes Lemon XL	200pk	MD IIA	VAH/DGHM, EN 1040, EN 13727	VAH/DGHM, EN 1275, EN 13624	EN 14348 (high bio load)	BGA(RKI)/DVV	x	x	x	x	Adeno Viruses (EN 14476)
9006410	Eurosept Max Wipes XL non-alc.	200pcs	MD IIA	VAH/DGHM,	VAH/DGHM, EN 1275, EN 1040, EN 1276, EN 13697	EN 14348 (high bio load)	BGA(RKI)/DVV	x	x	x	x	
9006411	Eurosept Max Wipes XL non-alc.	200pcs	MD IIA	VAH/DGHM,	VAH/DGHM, EN 1275, EN 1040, EN 1276, EN 13697	EN 14348 (high bio load)	BGA(RKI)/DVV	x	x	x	x	
9006138	Maxima non-alc. disinfection Wipes	75pcs	MD IIA	VAH/DGHM, EN 1040, EN 1276, EN 13697	VAH/DGHM, EN 1275, EN 1650, EN 13697	—	BGA(RKI)/DVV	x	x	x	x	“Noro (EN 14476)”
9006139	Maxima non-alc. disinfection Wipes	75pcs	MD IIA	VAH/DGHM, EN 1040, EN 1276, EN 13697	VAH/DGHM, EN 1275, EN 1650, EN 13697	—	BGA(RKI)/DVV	x	x	x	x	“Noro (EN 14476)”
9006281	Maxima non-alc. Surface	5l	MD IIA	VAH/DGHM, EN 1040, EN 1276, EN 13697	VAH/DGHM,EN 1275, EN 1650, EN 13697	—	BGA(RKI)/DVV	x	x	x	x	“Noro (EN 14476)”
9006137	Maxima non-alc. Surface disinfection	750ml	MD IIA	VAH/DGHM, EN 1040, EN 1276, EN 13697	VAH/DGHM,EN 1275, EN 1650, EN 13697	—	BGA(RKI)/DVV	x	x	x	x	“Noro (EN 14476)”

Itemcode	Description / Applicable Area	Contents	Classification Medical Device Class / Biocide / Cosmetical / Technical product	Bacteria in accordance to	Candida albicans in accordance to	TbB in accordance to	Membrane covered viruses i.e. HBV, HIV, HCV) in accordance to	Rotavirus	MRSA	H1N1	H1N5	Others
9006133	Eurosept Max Handdisinfection gel	1l	Biocide	VAH/DGHM, EN 1500, prEN 12054, EN 12791	VAH/DGHM, EN 1275	EN 14348 (high bio load)	BGA(RKI)/DVV	—	x	x	x	
9006132	Eurosept Max Handdisinfection gel	500ml	Biocide	VAH/DGHM, EN 1500, prEN 12054, EN 12791	VAH/DGHM, EN 1275	EN 14348 (high bio load)	BGA(RKI)/DVV	—	x	x	x	
9006135	Eurosept Max Handdisinfection gel	5l	Biocide	VAH/DGHM, EN 1500, prEN 12054, EN 12791	VAH/DGHM, EN 1275	EN 14348 (high bio load)	BGA(RKI)/DVV	—	x	x	x	
9006131	Eurosept Max Handdisinfection gel	75ml	Biocide	VAH/DGHM, EN 1500, prEN 12054, EN 12791	VAH/DGHM, EN 1275	EN 14348 (high bio load)	BGA(RKI)/DVV	—	x	x	x	
9006121	Eurosept Max Handdisinfection liquid	150ml	Biocide	VAH/DGHM, EN 1500, prEN 12054, EN 12791	VAH/DGHM, EN 1275	EN 14348 (high bio load)	BGA(RKI)/DVV	—	x	x	x	
9006128	Eurosept Max Handdisinfection liquid	1l	Biocide	VAH/DGHM, EN 1500, prEN 12054,	VAH/DGHM, EN 1275, EN 1275, EN 12791	EN 14348 (high bio load)	BGA(RKI)/DVV	—	x	x	x	
9006126	Eurosept Max Handdisinfection liquid	250ml	Biocide	VAH/DGHM, EN 1500, prEN 12054, EN 12791	VAH/DGHM, EN 1275	EN 14348 (high bio load)	BGA(RKI)/DVV	—	x	x	x	
9006129	Eurosept Max Handdisinfection liquid	500ml	Biocide	VAH/DGHM, EN 1500, prEN 12054,	VAH/DGHM, EN 1275, EN 1275, EN 12791	EN 14348 (high bio load)	BGA(RKI)/DVV	—	x	x	x	
9006130	Eurosept Max Handdisinfection liquid	5l	Biocide	VAH/DGHM, EN 1500, prEN 12054,	VAH/DGHM, EN 1275, EN 1275, EN 12791	EN 14348 (high bio load)	BGA(RKI)/DVV	—	x	x	x	
9002483	EuroSept Spittoon Cleaner	750ml	MD IIA	No effectivity report - cleaning product								
9006413	Maxima Nano-Spittoon Foam cleaner	400ml	MD IIA	No effectivity report - cleaning product								
9002091	Eurosept Max Foam Cleaner	1l	MD I	No effectivity report - cleaning product								
9001578	EuroSept Thermo Neutralizer	5l	MD I	No effectivity report - cleaning product								
9006140	EuroSept Thermo Cleanser	5l	MD I	No effectivity report - cleaning product								
9001577	EuroSept Thermo Powder	6kg	MD I	No effectivity report - cleaning product								
9001575	EuroSept Thermo Rinse	1l	MD I	No effectivity report - cleaning product								
9008395	Autoclave Cleaner	1000 ml	technical product	No effectivity report - cleaning product								
9002488	EuroSept Tray Cleaner Liquid concentrate	1l	technical product	No effectivity report - cleaning product								
9002487	EuroSept Tray Cleaner Liquid concentrate	5l	technical product	No effectivity report - cleaning product								
9006136	Eurosept Max Handwash	500ml	Cosmetical product	No effectivity report - cleaning product								
9002594	Eurosept Max Wash Lotion	1	Cosmetical product	No effectivity report - cleaning product								
9001317	EuroSept Max Wash lotion	1	Cosmetical product	No effectivity report - cleaning product								
9002608	Eurosept Max Wash Lotion Sensitive	1l	Cosmetical product	No effectivity report - cleaning product								
9002609	EuroSept Max Wash lotion Sensitive	5l	Cosmetical product	No effectivity report - cleaning product								
9006270	SafeSept Handdisinfection liquid	150ml	Pharmaceutical	VAH/DGHM, EN 1500, prEN 12054, EN 1275, EN 12791	AH/DGHM, EN 1275	EN 14348 (high bio load)	BGA(RKI)/DVV	—	x	x	x	
9006273	SafeSept Handdisinfection liquid	1l	Pharmaceutical	VAH/DGHM, EN 1500, prEN 12054, EN 1275, EN 12791	AH/DGHM, EN 1275	EN 14348 (high bio load)	BGA(RKI)/DVV	—	x	x	x	
9006271	SafeSept Handdisinfection liquid	250ml	Pharmaceutical	VAH/DGHM, EN 1500, prEN 12054, EN 1275, EN 12791	AH/DGHM, EN 1275	EN 14348 (high bio load)	BGA(RKI)/DVV	—	x	x	x	
9006272	SafeSept Handdisinfection liquid	500ml	Pharmaceutical	VAH/DGHM, EN 1500, prEN 12054, EN 1275, EN 12791	AH/DGHM, EN 1275	EN 14348 (high bio load)	BGA(RKI)/DVV	—	x	x	x	
9006274	SafeSept Handdisinfection liquid	5l	Pharmaceutical	VAH/DGHM, EN 1500, prEN 12054, EN 1275, EN 12791	AH/DGHM, EN 1275	EN 14348 (high bio load)	BGA(RKI)/DVV	—	x	x	x	

**itemcode According Norms:**

9001288	EN1040 / EN1275 / EN14476
9001289	EN1040 / EN1275 / EN14476
9001294	EN1040 / EN1275
9001295	EN1040 / EN1275
9002098	no survey available
9001296	no survey available
9003331	EN13727 / EN14561 / EN13624 / EN14562 / EN14348 / EN14563
9003332	EN13727 / EN14561 / EN13624 / EN14562 / EN14348 / EN14563
9002600	EN13727 / EN14561 / EN13624 / EN14562 / EN14348 / EN14563
9002601	EN13727 / EN14561 / EN13624 / EN14562 / EN14348 / EN14563
9001291	EN14348 / EN1040 / EN1275 / EN13624 / EN13727
9001290	EN14348 / EN1040 / EN1275 / EN13624 / EN13727
9002092	EN14348
9002089	EN14348 / EN1040 / EN1275 / EN13624 / EN13727
9002090	EN14348 / EN1040 / EN1275 / EN13624 / EN13727
9001293	EN14348 / EN1040 / EN1275 / EN13624 / EN13727
9001292	EN14348 / EN1040 / EN1275 / EN13624 / EN13727
9002095	EN14348 / EN1040 / EN1275 / EN13624 / EN13727
9002096	EN14348 / EN1040 / EN1275 / EN13624 / EN13727
9006410	DIN ISO 9001:2001 / DIN EN ISO 13485:2003; 93/42/EWG; DIN 12791 European Pharmaceutical book 2.2.6; DIN 1262; „Standardmethods of DGHM for checking and approval of chemical disinfectionmethods“ (Stand: 01.09.2001); Testmethod according BGA/DVV (1982); DIN EN ISO 2871-2; DIN EN 14971: 2006; 98/8/EG; EN 1041:1998”
9006411	“DIN ISO 9001:2001 / DIN EN ISO 13485:2003; 93/42/EWG; DIN 12791; European Pharmaceutical book 2.2.6; DIN 1262; „Standardmethods of DGHM for checking and approval of chemical disinfectionmethods“ (Stand: 01.09.2001); Testmethod according BGA/DVV (1982); DIN EN ISO 2871-2; DIN EN 14971:2006; 98/8/EG; EN 1041:1998”
9006138	“DIN ISO 9001:2001 / DIN EN ISO 13485:2003; 93/42/EWG; DIN 12791; European Pharmaceutical book 2.2.6; DIN 1262; „Standardmethods of DGHM for checking and approval of chemical disinfectionmethods“ (Stand: 01.09.2001); Testmethod according BGA/DVV (1982); DIN EN ISO 2871-2; DIN EN 14971:2006; 98/8/EG; EN 1041:1998”
9006139	“DIN ISO 9001:2001 / DIN EN ISO 13485:2003; 93/42/EWG; DIN 12791; European Pharmaceutical book 2.2.6; DIN 1262; „Standardmethods of DGHM for checking and approval of chemical disinfectionmethods“ (Stand: 01.09.2001); Testmethod according BGA/DVV (1982); DIN EN ISO 2871-2; DIN EN 14971:2006; 98/8/EG; EN 1041:1998”
9006281	“DIN ISO 9001:2001 / DIN EN ISO 13485:2003; 93/42/EWG; DIN 12791; European Pharmaceutical book 2.2.6; DIN 1262; „Standardmethods of DGHM for checking and approval of chemical disinfectionmethods“ (Stand: 01.09.2001); Testmethod according BGA/DVV (1982); DIN EN ISO 2871-2; DIN EN 14971:2006; 98/8/EG; EN 1041:1998”
9006137	“DIN ISO 9001:2001 / DIN EN ISO 13485:2003; 93/42/EWG; DIN 12791; European Pharmaceutical book 2.2.6; DIN 1262; „Standardmethods of DGHM for checking and approval of chemical disinfectionmethods“ (Stand: 01.09.2001); Testmethod according BGA/DVV (1982); DIN EN ISO 2871-2; DIN EN 14971:2006; 98/8/EG; EN 1041:1998”

**itemcode According Norms:**

9006133

9006132

9006135

9006131

9006121

9006128

9006126

9006129

9006130

9002483 No used norm because it is a cleaning product

9006413 No used norm because it is a cleaning product

9002091 EN1040 / EN1275

9001578 no survey available

9006140 no survey available

9001577 no survey available

9001575 no survey available

9008395 No used norm because it is a cleaning product

9002488 No used norm because it is a cleaning product

9002487 No used norm because it is a cleaning product

9006136 No used norm because it is a cosmetical product

9002594 no survey available

9001317 no survey available

9002608 no survey available

9002609 no survey available

9006270

9006273

9006271

9006272

9006274