

New products

Barrier sleeves

Barrier film clear, 100 x 150 mm, box 1200 pcs	1009869
Light handle cover, 100 x 150 mm, box 500 pcs	1009864
Air water syringe barrier sleeve, 60 x 250 mm, box 500 pcs	1009860
Chair cover, 700 x 500 mm, box 225 pcs	1009863



Available May 2013

Barrier sleeve, fits 7- and 8-mm- light guide, 400 pcs / box	9009042
Handpiece-barrier sleeve, 38 mm x 203,2 mm (1.5" x 8"), 500 pcs / box	9009043
Barrier sleeve for halogen light units, fits CU 1000, 250 pcs / box	9009044
Barrier sleeve for LED light units, 250 pcs / box	9009045
Barrier sleeve for tubes, 10,16 cm, 365 m / roll	9009046



Helixtest

Kit contains of 1 test body and 250 test strips	9008649
Refill: 250 test strips	9008650
Refill: 1 test body	9008651

Helix Test for steam sterilization validates every sterilization process in a loaded autoclave that air and other non-condensable gases have been removed. The indicator ink is non-toxic and the adhesive does not contain latex. Complies to EN 867-5, ISO 11140-1, Class 2

Disposable Propy Angle - Acclean® Latex Free

Webbed cup, soft cup, green, box with 100 pcs
Webbed cup, hard cup, white, box with 100 pcs
Disposable propy angles, ergonomically designed for ease of use.
Webbed cup, latex free. Individually wrapped,
100 pieces per box. Best to use with a straight handpiece such as BA 170 Ultimate

9008926
9008927





Handtowel dispenser for C-Fold & Z-Fold & hand roll towels

Made of high value resin. Colour: pearl white

Hand towel dispenser, big, Sizes: 374x275x140 9008446 for 9003324/9003325/9003326/9003327
Hand towel dispenser, small, Sizes: 374x275x106 9008447 for 9003324/9003325/9003326/9003327

Hand towel dispenser, big, Sizes: 365x250x230 9008448 for hand towels 9003315 / 9003316
Hand towel dispenser, small, Sizes: 345x180x172 9008449 for hand towels 9003314 / 9003318



Hot- and cold compresses

Sizes: 14 x 13 cm, box with 10 pcs 9008294
Sizes: 29 x 12 cm, box with 10 pcs 9008295
Sizes: 26,5 x 13 cm, box with 10 pcs 9008296

Re-usable. Provides effective relief from muscular pain and helps reduce swelling.
Great for treating a wide range of soft tissue injuries, pack will remain flexible at low temperatures. Pack lasts longer and is more effective than a standard ice or hot pack.
Reusable pack is perfect for continuous treatment at home. Cold: up to -20°C, hot up to +80°C.





New products

Replacement silicone rails for instrument cassette

bag with 3 pcs

5	Blue	9008829
	Red	9008830
	Green	9008831
	Violet	9008832
	Yellow	9008833
10	Blue	9008834
	Red	9008835
	Green	9008836
	Violet	9008837
	Yellow	9008838
20	Blue	9008839
	Red	9008840
	Green	9008841
	Violet	9008842
	Yellow	9008843



Gelantine Sponges

box with 50 pcs 9008806


Gelatine sponges with colloidal silver are an efficient and biocompatible product for post extraction treatment. The colloidal silver also has a bactericidal depot effect. The gelatine sponges promote coagulation and haemostasis, support healing, prevent infection through their antimicrobial effect, prevent post operative haemorrhage and the formation of secondary cavities. The sponges are resorbable, sterile and can be cut to size. Composition: One gelatine sponge (14 x 7 x 7 mm) contains: Hardened gelatine Ph. Eur. 9.5 mg Colloid silver Ph. Eur. 0.5 mg



Disposable filter for suction device

yellow, box 12 pcs

9006472



Phosphor plates envelopes

Size 0, 22 x 35 mm, box with 100 pcs 9008443
Size 1, 24 x 40 mm, box with 100 pcs 9008444
size 2, 31 x 41 mm, box with 300 pcs 9008445

Henry Schein Phosphor plates envelopes protect phosphor plate from exposure to light. Works with all brands of phosphor plate readers.
Does NOT require scissors to open up. Easy tear-away seam. Easy peel glue strip cover



Two additional colours:

Pink 9008845
Violet 9008846

Product details can be easily downloaded

Detailed product information on www.henryscheinbrand.com

latest product flyer on <http://www.henryscheinbrand.com/index.php?page=news>

Additional products

Prophy cups - latexfree - box 100 pcs
blue, hard: 9008907, Violet, medium: 9008908
Pink, soft: 9008909
Rubber Cups, screw type

Prophy brushes screwtype, latex free - box 100 pcs
Blue, hard: 9008904, Violet, medium: 9008905, Pink, soft: 9008906
Screwtype

Natural / Nylon Bristle Brush - box 100 pcs
Nylon brushes 9008911
Natur brushes 9008912
Screwtype

Product changes

HS-Bleaching

Products of the US American product line have been discontinued! Home Bleaching is now a cosmetic product
9008873 Bleaching home kit 10% CP (= 3.6%HP)
9008874 Bleaching home kit 16% CP (= 5.7% HP)
We discontinued all 22% bleaching products due to the new regulations. The In-Office-Bleaching-product stays medical device (new item code: 9008875). We will sell off the current stock levels of the remaining in office bleaching.



Scalpels and scalpel blades:

Scalpels and scalpel blades will be changed. In the future we are only going to have carbon steel products. Products have been tested by the internal dentist test panel.

Rapid Cold Pack:
old 9003330 — new 9008442



Desensitizing mousse

Kit: 2 syringes with 1,2 ml and 5 application tips 9008817

Henry Schein Desensitising Mousse provides intensive care for teeth also following a tooth whitening procedure. It remineralises the teeth and can relieve possible sensitivity. The cream has a pleasant fruity taste. Henry Schein Desensitising Mousse contains Sodium fluoride, Potassium nitrate and 30% Xylitol. The white paste can be used both in whitening trays, and directly on the teeth. This product is intended for use either by a dental professional in the office or provided to the patient for home treatment of dentin sensitivity.

 HENRY SCHEIN®

N
V
T
2
0
1
3
Germany



- High quality products from selected manufacturers
- Best value for your money
- Unconditional Guarantee

Henry Schein Own Brand
Quality products from a company you trust
www.henryscheinbrand.com

NVT 2013